
## **T2 Basic Tray Long**

## T2 Metal Trays: Content and Tray Layout



**Osteosynthesi** 



## **T2 Basic Tray Long**

## T2 Metal Trays: Set Configuration and Inventory Control Shee

	REF	Description	Quantity	Inventory	Notes
	T2 Basic Long				
-	702429	Elastosil Teardrop Handle, AO Coupling	1		
7-	703165	Protection Sleeve, Retrograde	1		
	1806-0022	Guide Wire Ruler	1		
5	1806-0032	Awl Plug	1		
	1806-0041	Awl, Curved 9mm **	1		
	1806-0110	Universal Rod, 9mm	1		
	1806-0125	Reduction Spoon, 9mm	1		
	1806-0130	Wrench 8mm/10mm	1		
	1806-0135	Insertion Wrench, 10mm	1		
1	1806-0150	Strike Plate	1		
	1806-0170	Slotted Hammer	1		
	1806-0170	Tissue Protection Sleeve, Long	1		
			1		
	1806-0203	Screwdriver, Self-Holding, Extra Short (hex 3.5)	1		
	1806-0215	Drill Sleeve, Long	1		
	1806-0227	Screwdriver Shaft AO, Long (hex 3.5)	<u>l</u>		
	1806-0233	Screwdriver, Self-Holding, Long (hex 3.5)	1		
	1806-0268	Screwdriver Shaft, Compression (hex 3.5)	1		
-	1806-0271	Guide Wire Pusher	1		
	1806-0315	Trocar, Long	1		
	1806-0325	Screw Gauge, Long (20 - 120mm)	1		
	1806-0331	Screw Gauge (20 - 120mm)	1		
	1806-0350	Extraction Rod, Conical	1	<u> </u>	
	1806-0365	Screw Scale, Long	1		
	1806-1095	Guide Wire Handle	1		
=	1806-1096	Guide Wire Handle Chuck (2-35mm)	1		
	1806-2014	Rigid Reamer Ø12mm	1		
	1806-9900	T2 Basic Long Instrument Tray (w. Lid)	1		
	Kit No. 1806-9901	<sup>2</sup> Basic Long Instrument Tray weighs 7.9 kg. is delivered fully equipped with instruments according to the l able in all countries.	list above.		
		nated "Outside of the U.S." may not be ordered for the U.S. ma	rket.		
		nation about the function and use of the T2 systems please refe erative Technique Guides.			

\*\* Either a 1806-0041 Ø9mm, or 1806-0040 Ø10mm, Awl Curved or 1806-0045 Awl, Straight, Ø10mm can be stored in the tray.