


T2 Basic Tray Short

T2 Metal Trays: Content and Tray Layout




Osteosynthesis



T2 Basic Tray Short

T2 Metal Trays: Set Configuration and Inventory Control Sheet

	REF	Description	Quantity	Inventory	Notes
	T2 Basic Short				
	702429	Teardrop Handle, AO Coupling	1	<input type="checkbox"/>	
	1806-0022	Guide Wire Ruler	1	<input type="checkbox"/>	
	1806-0032	Awl Plug	1	<input type="checkbox"/>	
	1806-0041	Awl, Curved 9mm **	1	<input type="checkbox"/>	
	1806-0113	Universal Rod, Short	1	<input type="checkbox"/>	
	1806-0130	Wrench 8mm/10mm	1	<input type="checkbox"/>	
	1806-0135	Insertion Wrench, 10mm	1	<input type="checkbox"/>	
	1806-0150	Strike Plate	1	<input type="checkbox"/>	
	1806-0170	Slotted Hammer	1	<input type="checkbox"/>	
	1806-0180	Tissue Protection Sleeve, Short	1	<input type="checkbox"/>	
	1806-0203	Screwdriver, Self-Holding, Extra Short (hex 3.5)	1	<input type="checkbox"/>	
	1806-0210	Drill Sleeve, Short	1	<input type="checkbox"/>	
	1806-0222	Screwdriver Shaft AO, Short (hex 3.5)	1	<input type="checkbox"/>	
	1806-0238	Screwdriver, Self-Holding, Short (hex 3.5)	1	<input type="checkbox"/>	
	1806-0271	Guide Wire Pusher	1	<input type="checkbox"/>	
	1806-0310	Trocar, Short	1	<input type="checkbox"/>	
	1806-0330	Screw Gauge, Short	1	<input type="checkbox"/>	
	1806-0353	Extraction Rod, Conical	1	<input type="checkbox"/>	
	1806-0360	Screw Scale, Short	1	<input type="checkbox"/>	
	1806-0363	Reduction Rod	1	<input type="checkbox"/>	
	1806-0390	Depth Gauge (Standard Style for Freehand Locking, 20 - 60mm)	1	<input type="checkbox"/>	
	1806-0410	Rigid Reamer Sleeve Ø10mm	1	<input type="checkbox"/>	
	1806-0411	Rigid Reamer Trocar Ø10mm	1	<input type="checkbox"/>	
	1806-1095	Guide Wire Handle (2 - 35mm)	1	<input type="checkbox"/>	
	1806-1096	Guide Wire Handle Chuck	1	<input type="checkbox"/>	
	1806-2020	Crown Drill	1	<input type="checkbox"/>	
	1806-9905	T2 Basic Short Instrument Tray (w. Lid)	1	<input type="checkbox"/>	
<div>Remarks:</div> <div>A fully equipped MMT T2 Basic Tray Short incl. Instruments and Drill Rack weighs 7.8 kg.</div> <div>Kit No. 1806-9906 is delivered fully equipped with instruments according to the list above –without Drill Rack. Listed Drills and K-Wires are delivered in sterile packaging. It may not be available in all countries.</div> <div>Instruments designated “Outside of the U.S.” may not be ordered for the U.S. market.</div> <div>For detailed explanation about the function and use of the T2 systems please refer to the corresponding Operative Technique Guides.</div> <div>** Either a 1806-0041 Ø9mm, or 1806-0040 Ø10mm, Awl Curved or 1806-0045 Awl, Straight, Ø10mm can be stored in the tray</div>					
<div>A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.</div> <div>The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.</div> <div>Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, T2. All other trademarks are trademarks of their respective owners or holders.</div> <div>Literature Number: B1000076 Rev. 0_08/10</div> <div>Copyright © 2010 Stryker</div> <div>Stryker Trauma GmbH Prof.-Küntscher-Straße 1 - 5 D - 24232 Schönkirchen Germany</div>					